

GeneFect

BROAD SPECTRUM DISINFECTANT



GEN100 - 1 Gallon 4/cs
Closed Loop

**Broad Spectrum
Quaternary Ammonium
Disinfectant/Cleaner**

Kills SARS-CoV-2 in 2 Minutes*

HMIS Health Rating - 2 (In-use)

Alkaline pH

2 ounce/Gallon Dilution

Fresh Scent



**Virus causing CoVid-19 disease*

Genesan

www.cleaneasier.com






GeneFect: Disinfectant/Cleaner



GEN100

GEN217

Product and Ordering Information

Item #	Packaging	Y	N	Case Size/Gallons	Dilution Ratio
GEN100	Closed Loop <input checked="" type="checkbox"/> Tipped <input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	 <p>4 Gal</p>	 2 oz/gal. (1:64)  2 oz/gal. (1:64)
pH (conc): 9.0 pH (1:64): 9.0 Fragrance: fresh Color: green Dilution: Extreme Units				Yield/Case: 260 Gal 260 Gal	
GEN217	Closed Loop <input checked="" type="checkbox"/> Tipped (1:64) <input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	 <p>2 Gal</p>	 2 oz/gal. (1:64) GEN217 not recommended for bucket-filling
Dilution: Extreme Units PODIUM.4 White Trigger (GEN237)				Yield/Case: 130 Gal	